

# CAGE-TB



## Automated smartphone-based cough audio classification for rapid tuberculosis triage testing (cough audio triage for tb)


The CAGE-TB app aims to modify TB diagnosis by using algorithms to differentiate between TB and non-TB coughs, identifying TB susceptible individuals and prioritizing them for additional sputum testing. The app would allow for efficient testing, earlier TB detection, and reducing transmission.



For more information:



The Cough Audio Triage for Tuberculosis (CAGE-TB) app. a person will need to cough into their smartphone’s speaker. They will then receive their diagnosis, whereafter they will receive specialised care. INFOGRAPHIC: Marianne Francis Stewart, 2021.



**Figure 1: Cage-TB app.** Initial screen of the app presents four options. Tapping ‘start now’ directs users to a second screen, where they can input the patients’ ID for a new record. A red recording indicator signifies whether a recording is in progress.



For more information on AI and robotics in healthcare:

